


YOUR SMILE
IS OUR
EVERYTHING



K-A II High

Clinic Line > Dental Amalgam > Capsule Amalgam

High silver alloy which contains a mixture of spherical and irregular particles. Predosed in self activating capsules. Gamma-2 free. It has 3 spill options as spill 1, spill 2 and spill 3

And also powder Amalgam %45 & %69

Indications

Stress bearing fillings in posterior teeth (Class I and II)

Retrograde root canal filling after apex resection

Features

Mixture of spherical and irregular particles

High silver content

Exceptionally durable

Easy to process

Relatively fast hardening.



Composition

1- %045 Ag %30,5 Sn %24 Cu 5,04 Zn

2- %55 Ag %29 Sn %15,8 Cu %0,3 Zn

3- %69,2 Ag %18,6 Sn %11,9 Cu %0,3 Zn

Packaging

50 or 500 capsule in a jar

50 capsule in a carton box

50 or 100 capsule in a nylon bag

30 grams powder amalgam in a jar



K-Compo

Clinic Line > Dental Composite > Nano Hybrid Composite

K-Compo is a radio opaquemicrohybride light cure composite system for anterior and posteriorfield (Extra white & Posterior)

Properties

Universal Composite

It contains resins based Bis GMA

It has high filler rate

It has good polishing facility

It has high resistance to abrasion

It provides ease of use

It is radio opaque

It provides ease of use

Non stick to instruments

Packaging

Single

4,0 gr syringe


Kit

4,0 gr syringe x 12 or 10 or 8 or 6 pcs

5 ml K-Sole Bond single bonding x 1 pcs

3 ml K-Etch x 1 pcs. and Accessoires





K-Advantage

Clinic Line > Dental Composite > Nano Hybrid Composite

K-Advantage is radio opaque nano composite system which can be photocurable, applied in the posterior and anterior field (A1-A2-A3-A3.5-B1-B2-B3-C1)

Properties	Low shrinkage
It contains Bis- GMA, UDMA, Bis- EMA resins	Non sticky
It has high filler rate	Packaging
It can be polished with high degree	Single
It has high resistance to abrasion	4,0 grsyringe
It is radio opaque	Kit
It provides ease of use	4,0 grsyringe x 12 or 10 or 8 or 6 pcs.
Long durability	5 ml K-Sole Bond single bonding x 1
Higher hardness	5 ml K-SEBond self-etch bonding x 1 pcs.



K-Advantage Nano

Clinic Line > Dental Composite > Nano Hybrid Composite

K-Advantage Nano is radio opaque nano composite system which can be the photocurable, applied in the posterior and anteriorfield (A1-A2-A3-A3.5-B1-B2-B3-C1)

Properties	Low shrinkage
It contains Bis- GMA, UDMA, Bis- EMA resines	Non sticky
It has high filler rate	Packaging
It can be polished with high degree	Single
It has high resistance to abrasion	4,0 grsyringe
It is radio opaque	Kit
It provides ease of use	4,0 grsyringe x 12 or 10 or 8 or 6 pcs.
Long durability	5 ml K-Sole Bond single bonding x 1
Higher hardness	5 ml K-SEBond self-etch bonding x 1 pcs.
Non sticky	3 ml K-Etch x 1 pcs. and Accessoires



K-Etch

Clinic Line > Bonding Etching

K-Etch is an etching gel containing %37 phosphoric acid. Delivered in 3,0 ml.syringe with very thin application needles.

Features	Easily Rinse
Blue coloured	Packaging
Perfect Surface Retention	3 ml syringe x 2



Coming soon...

K-Polishing Disk & K-Polishing Strips

Clinic Line > DentalPolishing Disk



K-matrix

CLINIC LINE > DENTALMATRIX

It is the best way to make smart and simple posterior applications. It has special design with transparent or steel and two different type of height. It is important to keep patient comfort during treatment.

Packaging

50 pcs. n a box

Size

- Metal adapt Smart Matris thickness 0,03 height 5.0
- Metal adapt Smart Matris thickness 0,03 height 6,3
- Metal adapt Smart Matris thickness 0,05 height 5.0
- Metal adapt Smart Matris thickness 0,05 height 6.3
- Metal adapt Smart Matris thickness 0,038 height 5.0
- Metal adapt Smart Matris thickness 0,038 height 6,3





K-Sole Bond

Clinic Line > Bonding Etching

K-Sole Bond is light curing, single component bonding agent for the adhesive restorations. It is universal and compatible with all composites.

Properties

Properties

Elastic connection

Packaging

5 ml in a bottle



K-SE Bond

Clinic Line > Bonding Etching

K-SE Bond is light curing, single component bonding agent for the adhesive restorations. It is universal and compatible with all composites.

Properties

Properties

Elastic connection

Packaging

5 ml in a bottle






Rotary File

Clinic Line > Dental Composite > Nano Hybrid Composite

Endo-Art launched %100 Made In Turkey Endodontic Rotary File system. Heat treated control memory nickel-titanium file system. It is perfectly designed for great endodontic treatments.

Properties	Non stick to instruments
Excellenot Flexibility	Packaging
Perfect Strength	6 file system in blister
Non-cutting tip	
21 mm, 23 mm, 25 mm, 27 mm and 31 mm	



EndoArt H & K Files Niti

Clinic Line > DentalComposite

Endodontic hand use H type and K type files with excellent cutting efficiency. The high quality raw material is Nickel Titanium .

Properties	Excellent handling
It has silicon stops	Packaging
Identifically Numbering	10 pcs. blister in a box
Codes: #06/08/10/15/20/25/30/35/40/45/50/55/60/70/80/90/100/110/120/130/140/	
onforms IS 33 -1 standarts	
Sterilization 134 C	
21 mm, 25 mm, 31 mm length	
Ergonomic design	



+90 242 324 9092



+90 546 454 9745



info@kuantadental.com



kuanta.group@gmail.com



www.kuantadental.com



Zümrütova Mah. Sinanoğlu Cad.
Osman Eceoğlu Apt. No:33.1 / Muratpaşa / ANTALYA

